






WORKING AT A CELLULAR LEVEL TO ENHANCE ALL ASPECTS OF OPTIMAL HEALTH & HEALING

F R E Q U E N C Y S E T T I N G S

PROGRAM	NAME OF SETTING	DESCRIPTION OF SIGNALS Basic Signal	User Signal	DEFAULT TIME	DEFAULT INTENSITY	SLOPE OF INTENSITY
1	 Increase Heart and Blood Circulation	100%	Beta to Alpha (downwards)	15 minute treatment	1	↓
2	 Relaxation	10%, decreases to 5%	Beta to Alpha to Theta, to Delta (downwards)	15 minute treatment	3	↓
3	 Energy	100% to Beta	Theta to Beta (upwards)	15 minute treatment	5	↑
4	 Decrease Inflammation/ Pain Relief	100% downward in steps of 2HZ	Beta to Theta (downwards)	15 minute treatment	3	↓
					Strongest Program Use Default Settings Only!	
5	 Recovery/Prevention	15 minutes Program 2, 15 minutes Program 3	Beta to Alpha to Theta, to Delta (downwards 2-1/2 min.) alternating w/Theta to Beta (upwards 2-1/2 min.)	30 minute treatment	6	↓ Then ↑

User's Signal Frequencies Affecting Electromagnetic-Cell-Communication: Directs the body of the horse.

Beta: 13-32 Hz - Alertness / Alpha 8-13 Hz - Relaxed / Theta: 5-8 Hz - Refresh Brain / Delta: 1-5 Hz Regenerate Brain

Basis Signal Frequencies Affecting Cellular Metabolism: Improved calcium effect and ion transfer within the cells.

Basis Signal Frequencies (300 - 600 - 800 - 1000 Hz)